

# AsTeRICS Workshop – HandsOn Part

---



# Welcome !!

- In this Workshop we will use the AsTeRICS System
- First we make sure that everybody has ARE and ACS installed
- Then we will build some easy models together
- After that you can choose an Assistive Technology project and build it together in small groups

# The Basics

---



- AsTeRICS can be downloaded from the official site:  
<http://www.asterics.org>
- The latest source code is available at GitHub  
<https://github.com/asterics/AsTeRICS>
- ACS is the graphical editor, which sends models with connected plugins to the ARE (the runtime environment)
- Documentation is available:
  - User Manual with step-by-step model creation guide
  - Developer manual for creating new plugins
  - Use **F1** in ACS when plugin is selected to display plugin-help !
- Let's start ACS and ARE and have a look !

# Some useful plugins – short examples:

---

- Sensors

- ComputerVision: FaceTrackerLK
- Generic: **DigitalIn**, **AnalogIn**
- GUI: ButtonGrid, Slider, Cellboard
- InertialMeasurement: **Acceleration**
- Simulation: SignalGenerator, Timer

- Processors

- Audio&Voice: **Speechprocessor**
- BasicMath: MathEvaluator, Threshold, Averager, Constant Dispatcher, Comparator, Differentiate / Integrate
- DataConverters: IntToString
- SignalShaping: Deadzone, SignalTranslate
- Event&StringProcessing: StringDispatcher
- Microcontroller: **Arduino**
- OSKA
- SignalPathways: PathSelector

- Actuators


- Audio: WavePlayer, SpeechSynthesizer
- GenericOutput: **DigitalOut**
- GUI: Oscilloscope, BarDisplay, EventVisualizer
- HomeControl: **FS20**, **IRTrans**
- InputDeviceEmulation: Mouse, Keyboard

\* Needs dedicated HW or SW / driver

# Example 1: Mouse-Control by Head Movements

- In this example we will show how you can control the mouse cursor via head movements
- We use a webcam and the FacetrackerLK sensor plugin (which tracks face movements)
- Mouse clicks are generated via a period of inactivity (no movements)

Note: you find a similar description in the AsteRICS User manual



just  
another  
example

# Example 1: Mouse-Control by Head Movements

---

## How to provide headtracking-controlled mouse alternatives for computer input

- The x- and y-position of the local mouse will be controlled by the user's head movement

### Requirements:

- A webcam


### Remarks:

- Camera position directly in front of the user, distance 70-120cm
- No other persons face should be in the field-of-view of the camera




## Example 1: Mouse-Control by Head Movements

- Create a new model in the ACS (System → **New Model**)
- Insert the „**FacetrackerLK**“ sensor component  
(Components → Sensors → ComputerVision → FacetrackerLK)
- Insert the „**Mouse**“ actuator  
(Components → Actuators → Input Device Emulation → Mouse)
- Connect ports „noseX“ to „mouseX“ and „noseY“ to „mouseY“



# Example 1: Mouse-Control by Head Movements

- Mouse plugin properties:
  - adjust **xMax** and **yMax** to the desired values, e.g. screen resolution
  - enter 0 for xMax and yMax to enable automatic detection of screen resolution
  - **deselect** „absolutePosition“  
this defines that the X and Y input values are relative changes → fits the output of the Facetracker component



Properties (Ctrl-P)

**Component:** Mouse.1

Search

Internal Properties


Component Name	Mouse.1
Component Type	asterics.Mouse
Component Descriptio	Local mouse ac
Component Class	actuator

Properties

enableMouse	<input checked="" type="checkbox"/>
absolutePosition	<input type="checkbox"/>
xMin	0
xMax	1366
yMin	0
yMax	768

# Example 1: Mouse-Control by Head Movements


- Mouse plugin, Input Port rider:



- Select „**Synchronize Inputs**“ option in the Input Port Riders for both inputs (mouseX and mouseY)

→ this will wait for both input coordinates to arrive before the mouse position is updated


- You can try out your model by uploading it to the ARE and clicking „Start Model“






## Example 1: Mouse-Control by Head Movements

- No clicking function implemented yet
- → further development: Add **dwel clicking**
- Adding processing component „**Deadzone**“  
(Components – Processors – Signal Shaping – Deadzone)
- Connect noseX/noseY outputs of the FacetrackerLK to the inX/inY inputs of the Deadzone




## Example 1: Mouse-Control by Head Movements

- Use the **Deadzone** component to define a desired movement level to start or stop the timing for the dwell click
- Deadzone component fade out x/y signal values in an adjustable range and generate event trigger if the x/y values are in- or outside this range
- Parameter „**radius**“ defines this range  
→ here it is the amount of nose movement
- Leave the radius at the default value of 10




# Example 1: Mouse-Control by Head Movements

- Default value 10 – movement range is set to 10 pixels from previous to current position
- Select „**Synchronize Inputs**“ option in the Input Port Riders of the Deadzone plugin for inX and inY



# Example 1: Mouse-Control by Head Movements


- How can we measure a certain timespan of low movement (to generate a dwell click) ?
- → Insert a **Timer** sensor component (Sensors – Simulation – Timer)
  - this component measures time, generates events if a time period has passed, perfoms timing loops
- Set time-period to 1000 in the components properties



# Example 1: Mouse-Control by Head Movements

- Connect event trigger port of the Deadzone component (purple) to the event listener of the Timer component (green)
- Click on the event channel (connection line is now purple dotted)

→ Events can be selected from dropdown menu




# Example 1: Mouse-Control by Head Movements

- Select „**enterZone**“ event for the „start“ function
- Select „**exitZone**“ event for the „stop“ and „reset“ function
  - These event connections control the Timer components
  - If nose movements stay below selected level of 10 pixels, the Timer is started
  - Else, the Timer is resetted to 0 and stopped
  - If the movement stays low for the full time period, the timer will generate its „**periodFinished**“ event.

Events (Ctrl-E)	
Timer.1	Deadzone.1
start	enterZone
start	---
stop	exitZone
stop	---
reset	exitZone
reset	---


# Example 1: Mouse-Control by Head Movements

- Draw a channel from the Timer's event trigger port (purple) to the event listener port of the Mouse (green)



## Example 1: Mouse-Control by Head Movements

- By clicking on the new channel (line is now purple dotted), the Events can be adjusted
- Assign the „periodFinished“ event to the „leftClick“ function
- After these settings, the model is usable and provides left-click


A screenshot of a software interface window titled "Events (Ctrl-E)". The window contains a table with two columns: "Mouse.1" and "Timer.1". The "Timer.1" column has a red oval around the "periodFinished" entry in the first row. The table lists various mouse events in the "Mouse.1" column and their corresponding functions in the "Timer.1" column.

Mouse.1	Timer.1
leftClick	periodFinished
leftClick	---
middleClick	---
rightClick	---
doubleClick	---
dragPress	---
dragRelease	---
wheelUp	---
wheelDown	---
center	---
activate	---
deactivate	---
toggle	---




# Example 1: Mouse-Control by Head Movements

**Include GUI elements for adjustable mouse acceleration and dwell timing as well as different click-actions**



# Example 1: Mouse-Control by Head Movements

- Add a **Slider** component (Components – Sensors – Graphical User Interface – Slider)
- Slider properties:
  - range of value can be defined (we can leave it at 0-100)
  - set slider's component name to "Mouse Speed"
  - set minorTickSpacing to "10"

A screenshot of a software development environment's Properties window. The window title is "Properties (Ctrl-P)". The component being edited is "Component: Mouse Speed". There is a search bar with a search icon and a search field containing the text "Search". The properties are organized into two sections: "Internal Properties" and "Properties".

Internal Properties	
Component Name	Mouse Speed
Component Type	asterics.Slider
Component Descriptio	a slider GUI element which provide
Component Class	sensor

Properties	
min	0
max	100
default	50
autosend	<input checked="" type="checkbox"/>
caption	Mouse Speed
majorTickSpacing	20
minorTickSpacing	10
alignment	horizontal
fontSize	14

## Example 1: Mouse-Control by Head Movements


---

- To modify the x/y mouse speed with the slider's values, a **MathEvaluator** processing component is needed (Components - Processors – Basic Math – MathEvaluator)
- First, the x-signal is modified:
  - Delete port connection from noseX to mouseX
  - Draw new port connection from “value” (Slider) to “inA” (MathEvaluator)
  - Draw new connection from noseX (Facetracker) to “inB” (MathEvaluator)
  - Draw a new connection from output port (MathEvaluator) to mouseX input port

# Example 1: Mouse-Control by Head Movements

- MathEvaluator properties:
  - Adjust “**expression**” property of the MathEvaluator – this defines what will be done with the inputs – in our case we will multiply inA and inB
  - Slider position <50 shall slow down mouse speed, slider positions >50 shall increase mouse speed

→  $a/50*b$




## Example 1: Mouse-Control by Head Movements

---

- For modifying the y-signal, copy and paste (Ctrl+C & Ctrl+V) the MathEvaluator
- New connections as for the x-direction:
  - Delete port connection from noseY to mouseY
  - Draw new connection from “value” (Slider) to “inA” of the second MathEvaluator
  - Draw new connection from noseY (Facetracker) to “inB” of the second MathEvaluator
  - Draw a new connection from output port of the second MathEvaluator to mouseY input port

# Example 1: Mouse-Control by Head Movements



# Example 1: Mouse-Control by Head Movements

---


## Adding different mouse click activities via GUI by

- adding a **ButtonGrid**
  - to select next click type
- informing mouse element about the next desired mouse click
  - by sending an “action string” to the Mouse element

**Action strings** contain commands which are understood by a number of specialized actuator elements. These string contain the addressed component and the desired command  
e.g. “@MOUSE:nextclick,right”

## Example 1: Mouse-Control by Head Movements

- Add the **ButtonGrid** component (Component – Sensors – Graphical User Interface – ButtonGrid)
- ButtonGrid properties:
  - Set “buttonCaption” properties of button 1, 2 and 3 to
    - RightClick
    - DoubleClick
    - DragClick
  - Choose “horizontalOrientation” for the ButtonGrid
  - Set a desired caption, e.g. “Next click is a”

A screenshot of a software development environment's Properties window for a ButtonGrid component. The window has a title bar that says "Properties" and a small arrow icon on the left. It contains a table with two columns: the first column lists properties and the second column shows their values. The properties listed are caption, horizontalOrientation, buttonCaption1, buttonCaption2, buttonCaption3, buttonCaption4, and buttonCaption5. The caption is "Next mouse click is a", horizontalOrientation is checked with an 'X' in a box, and the first three buttonCaption properties are set to "Right Click", "Double Click", and "Drag Click" respectively. The last two buttonCaption properties are empty.

Property	Value
caption	Next mouse click is a
horizontalOrientation	<input checked="" type="checkbox"/>
buttonCaption1	Right Click
buttonCaption2	Double Click
buttonCaption3	Drag Click
buttonCaption4	
buttonCaption5	



## Example 1: Mouse-Control by Head Movements

---

- Add the **StringDispatcher** component  
(Component – Processors – Event and String Processing – StringDispatcher)
  - Translates incoming events into outgoing strings
  - If buttons are pressed, desired actions strings are generated for the Mouse components
- Connect event trigger port of ButtonGrid (purple) to the event listener port of the StringDispatcher (green)
- Click on the event channel and attach
  - button1 to dispatchSlot1
  - button2 to dispatchSlot2
  - button3 to dispatchSlot3


## Example 1: Mouse-Control by Head Movements

---

- Define the strings for slot1-slot3 in the StringDispatcher properties:
  - slot1(button1): “@MOUSE:nextclick,right”
  - slot2(button2): “@MOUSE:nextclick,double”
  - slot3(button3): “@MOUSE:nextclick,drag”
- Connect output port of the StringDispatcher to the “action” input port of the mouse

## Example 1: Mouse-Control by Head Movements

- Change to GUI Designer tab to define a desired position for the ButtonGrid
  - Possible positions:




- "Upload Model" and "Start Model" to try it out!


# Example 1: Mouse-Control by Head Movements

## Further improvement tip


- If the mouse cursor is jittering a lot, use an averager (Components – Processors – Basic Math – Averager) with a BufferSize of 5 after Facetracker Outputs



# Example 1: Mouse-Control by Head Movements




## Example 2: Using the Arduino for interfacing




- Based on Atmel ATmega328
- 6 Analog Input Pins
- 14 Digital I/O Pins
- 32 KB Flash Memory
- 2 KB SRAM
- 1 KB EEPROM
- 16 MHz Clock Speed

- **Special Firmware** supports ARE – plugin
- Processors → Microcontroller Interface → Arduino



## Example 2: Using the Arduino for interfacing


- Installation of **COM Port**:
  - Attach Arduino to USB Port
  - Driver must be updated in Device manager to get COM Port
  - Select location of „Arduino.inf“ (folder: [CIMs/Arduino/driver](#))
  - Click „Install this driver software anyway“
- For Windows 8:
  - Reboot without Driver Signature:
    - Enter „shutdown /r /o /f /t 00“ in shell-window
    - Choose: troubleshooting -> advanced options
    - > startup settings -> reboot without driver signature enforcement



## Example 2: Using the Arduino for interfacing


---

- **Breadboard connections**
  - allow flexible building test circuits
  - the bars are connected vertically
  - the rails are connected horizontally
  - use cables to connect components and Arduino Pins







## Example 2: Connecting a LED to Arduino




- Connect the **Anode** of the Led to Pin7 (positive lead, usually the long leg)
- Connect the **Kathode** of the Led to a 470 Ohm resistor
- Connect the resistor's other end to GND (0 Volt)



## Example 2: Connecting a LED to Arduino


Microcontroller can apply  
5V (high) or 0V (low)  
to a digital output Pin



## Example 2: Connecting a LED to Arduino

- In the Arduino processor plugin's properties:
  - Define pin7Mode „output, default high“ or „output, default low“
  - Use Events Listeners to set output high or low (turn on/off)



Properties	
periodicADCUpdate	0
pin2Mode	not used
pin3Mode	not used
pin4Mode	not used
pin5Mode	not used
pin6Mode	not used
pin7Mode	output, default high
pin8Mode	not used
pin9Mode	not used
pin10Mode	not used
pin11Mode	not used
pin12Mode	not used
pin13Mode	not used



setPin6	---
clearPin6	---
setPin7	button1
setPin7	---
clearPin7	button2
clearPin7	---
setPin8	---
clearPin8	---
setPin9	---
clearPin9	---
setPin10	---


## Example 2: Connecting a PushButton to Arduino

- Connect one side of switch to GND and the other side to Pin 3



- When switch is pressed: all 4 leads connected !


## Example 2: Connecting a PushButton to Arduino



**Resistor pulls Pin-voltage to 5V (high)  
if the switch is not pressed**

**Button connects Pin to 0V (low)  
if pressed**


**Microcontroller can measure  
voltage (high or low)  
at a digital input Pin**



## Example 2: Connecting a PushButton to Arduino

- In the Arduino processor plugin:
  - Select pin3Mode: “input with pullup”
  - Event is triggered when PushButton is pressed (high -> low)
  - Event is triggered when PushButton is released (low -> high)



Properties	
periodicADCUpdate	0
pin2Mode	not used
pin3Mode	input with pullup
pin4Mode	not used
pin5Mode	not used
pin6Mode	not used
pin7Mode	not used
pin8Mode	not used
pin9Mode	not used
pin10Mode	not used
pin11Mode	not used
pin12Mode	not used
pin13Mode	not used



eventDisplay	pin3ChangedLowToHigh
eventDisplay	pin3ChangedHighToLow
eventDisplay	---


## Example 2: Switching a Led with a Pushbutton

- Switch is connected to Pin3 (internal PullUp)
- LED is connected to Pin7
- Switch Open → LED on
- Switch Closed → LED off



Properties	
periodicADCUpdate	0
pin2Mode	not used
pin3Mode	input with pullup
pin4Mode	not used
pin5Mode	not used
pin6Mode	not used
pin7Mode	output, default high
pin8Mode	not used
pin9Mode	not used
pin10Mode	not used
pin11Mode	not used
pin12Mode	not used
pin13Mode	not used

# Example 2: Switching a Led with a Pushbutton





Events (Ctrl-E) ▼ 🔍

Arduino.1	Arduino.1
setPin2	---
clearPin2	---
setPin3	---
clearPin3	---
setPin4	---
clearPin4	---
setPin5	---
clearPin5	---
setPin6	---
clearPin6	---
setPin7	pin3ChangedLowToHigh ▼
setPin7	---
clearPin7	pin3ChangedHighToLow ▼
clearPin7	---
setPin8	---




## Example 2: Reading Analogue Values

- Connect a Potentiometer
  - Right lead goes to 5 V
  - Middle lead goes to A0
  - Left leads goes to GND
- Voltage at A0 will vary from 0V to 5V when turning





## Example 2: Reading Analogue Values

- Arduino Analog-Digital Converter
  - 10 bit  $\rightarrow$  1024 Values
  - Zero Volt = Value 0
  - 5 Volt = value 1023
- Important: Set periodicADCUpdate to a positive Time (in milliseconds) otherwise the ADC is off !!!



Properties	
periodicADCUpdate	10
pin2Mode	not used
pin3Mode	not used
pin4Mode	not used
pin5Mode	not used
pin6Mode	not used
pin7Mode	not used
pin8Mode	not used
pin9Mode	not used
pin10Mode	not used
pin11Mode	not used
pin12Mode	not used
pin13Mode	not used

## Example 2: Reading Analogue Values




A voltage divider with a resistive sensor (e.g. LDR – Light dependent resistor) could also be a force sensor, temperature sensor etc.

Try: to print the brightness value in the ARE Window !


## Example 2: Reading Analogue Values

- Read a resistive sensor:
  - build a voltage divider!
  - R1 depends on the Sensor



## Example 2: Control Led with Sensor

- Try: Switch LED on if Sensor value reaches threshold




# Example 2: Control Led with Sensor

Component Description: Arduino microcontroller C++


Properties

periodicADCUpdate	50
pin2Mode	not used
pin3Mode	not used
pin4Mode	not used
pin5Mode	not used
pin6Mode	not used
pin7Mode	output, default low
pin8Mode	not used
pin9Mode	not used
pin10Mode	not used
pin11Mode	not used
pin12Mode	not used
pin13Mode	not used



Properties

condition	a greater threshold
outputMode	output a if condition met
eventMode	create event only if condition
threshold	800
threshold2	10



Events (Ctrl-E)

Arduino.1	Comparator.1
setPin2	---
clearPin2	---
setPin7	conditionTrue
setPin7	---
clearPin7	conditionFalse
clearPin7	---

# Hands-On Projects

---

**HANDS ON !!**



- Select a project and solve it in small groups !
- We have some sets of hardware which can be used
- 6 Project topics are available

# Project 1: Digitally controlled Mouse

---

- Realize a fully working mouse with 3 push buttons
  - 2 buttons: left/right or up/down  
pressing both buttons: selection of axis (X or Y)
  - 1 button: mouse click  
short press: left click  
long press: right click





# Project 2: Analog-controlled Mouse



---

- Realize a fully working mouse with 2 push buttons and 1 potentiometer
  - potentiometer: cursor movement with 2 directions (e.g. +x, -x)
  - 1 button: toggle axis (X or Y)
  - 1 button: mouse click
    - short press: left click
    - long press: right click



# Project 3: Virtual Keyboard

- Realize a keyboard with scanning via the Acceleration Sensor and/or EMG input
  - Use OSKA with a keyboard grid for writing
  - Use automatic scanning for key selection
  - Use EMG sensor for controlling the scanning
    - Minimal muscle movement selects key



# Project 4: Environmental Control

---

- Realize an environmental control system
  - Speech Recognition input
  - FS20 for controlling a 220V light
  - Abotic door opener connected via GPO module



# Project 5: Infrared Control

---


- Control a GhettoBlaster via IR device
  - GUI (Cellboard) with automatic scanning
  - IR Trans for infrared remote control
  - Send commands to HiFi Stereo Radio:
    - On/Off
    - Play / Stop
    - Volume control
    - Next/Previous song



# Project 6: Accessible Pong Game

---

- Realize a Pong Game Interface via desired input sensors
  - Player 1: uses Accelerometer: tilting controls paddle 1
  - Player 2: uses WebCam / FaceTracking: up/down movement of head controls paddle 2



# Project 1 – Tips


---

- Connect 3 buttons to your DigitalIn as shown:
- Alternatively, you can use an Arduino and use 3 pushbuttons



# Project 1 – Tips

- To generate 1 right-click with the Mouse actuator, you can send the string „@MOUSE:nextclick,right“ to the action input port
- Use a **StringDispatcher** processor to send desired strings to other plugins (see HeadMouse example)




# Project 1 – Tips

---

- Use a **Timer** sensor to detect a long button press
- If the button is held for a desired time (defined in the properties) the timer can send an event
- Attach Events from the DigitalIn to the Timer plugin so that:
  - Pressing the button starts the timer
  - Releasing the button stops (and resets) the timer.

The example shows events from a button connected to DigitalIn 3:


A screenshot of a software interface showing an event table. The table has two columns: "Timer.1.2" and "DigitalIn.1". The table contains six rows of event data. The first row shows a "start" event from the timer triggered by the "in3Low" event from the digital input. The second row shows a "start" event from the timer with no specific digital input event listed. The third row shows a "stop" event from the timer triggered by the "in3High" event from the digital input. The fourth row shows a "stop" event from the timer with no specific digital input event listed. The fifth row shows a "reset" event from the timer triggered by the "in3High" event from the digital input. The sixth row shows a "reset" event from the timer with no specific digital input event listed.

Timer.1.2	DigitalIn.1
start	in3Low
start	---
stop	in3High
stop	---
reset	in3High
reset	---



# Project 1 – Tips


- Use a ConstantDispatcher processor to create movement:
  - Send negative value (-1) to Mouse actuator x/y for up/left
  - Send positive value (1) to Mouse actuator x/y for right/down



# Project 1 – Tips


---

- Dispatch up/left or right/down movement via connected pushbuttons (here: DigitalIn Pins)

A screenshot of the 'Events (Ctrl-E)' window in a software IDE. The window is divided into two columns: 'ConstantDispatcher.1' and 'DigitalIn.1'. The 'DigitalIn.1' column contains a list of digital input pins with dropdown menus. The pins are: in5Low, ---, in6Low, ---, in5High, in6High, and ---. Each pin has a small downward arrow next to it.


ConstantDispatcher.1	DigitalIn.1
dispatchSlot1	in5Low ▼
dispatchSlot1	--- ▼
dispatchSlot2	in6Low ▼
dispatchSlot2	--- ▼
dispatchSlot3	in5High ▼
dispatchSlot3	in6High ▼
dispatchSlot3	--- ▼

# Project 1: a possible solution ...



# Project 2 – Tips


---



# Project 2 – Tips

---

- Different value ranges of mouse movement and Arduino ADC:
  - Mouse: relative, positive and negative values
  - ADC: absolute, only positive values (0-1023, 12bit)
- Solution: SignalTranslation plugin

A screenshot of a software interface showing the properties of a component named "SignalTranslation.1". The interface includes a search bar, a list of internal properties, and a list of general properties.

Component: SignalTranslation.1

Search

Internal Properties

Component Class	processor
Component Name	SignalTranslation.1
Component Type	asterics.SignalTranslation
Component Description	Translates signals from a certain input range to an output range


Properties

inMin	0
inMax	1023
outMin	63
outMax	-64


# Project 2 – Tips

---

- It's hard to find the exact centre of the potentiometer
  - > the Cursor is always moving
- We need a window around the centre, where the cursor does not move:
  - > Deadzone plugin



# Project 2: a possible solution ....



# Project 3 – Tips

---

- Connect the EMG shield to the Arduino
  - Attach electrodes to forearm
- Use an Arduino plugin to view analog values of AD0
  - Select an update rate of 5ms (200Hz) for the Arduino ADC
- Use oscilloscope to view raw signal (try it out !)



# Project 3 – Tips


---

- To calculate the EMG force:
  - With a **Differentiate** processor you to get rid of the DC-offsets and slow changes in the EMG signal. (You can also experiment with the Filter Plugin)
  - use a **MathEvaluator** to get the absolute values of the signal enter **abs(a)** as expression in the MathEvaluator properties
  - use an **Averager** to calculate the average of 30 samples
- Use a **Threshold** processor to create events if the EMG force exceeds a certain level
  - Check reasonable values with oscilloscope or BarDisplay

# Project 3 – Tips


---

- Add a „**OSKAInternalScanning**“ plugin
  - this launches the OSKA on-screen keyboard in automatic scanning mode
  - Upload model to ARE (to update list of available keyboards)
  - Select an appropriate keyboard grid for writing for example: „**keyboards\writing\abc.xml**“
  - Select scan speed 4 and highlightStyle „**InvertKeys**“

A screenshot of the 'Properties' dialog box for the OSKA Internal Scanning plugin. The dialog has a title bar with a small triangle icon and the text 'Properties'. It contains a list of properties with their corresponding values and controls.

Property	Value	Control
port	4546	Spin box
title	OSKA keyboard	Text field
oskaPath	..\oska\start keyboard.exe	Text field
keyboardPath	keyboards\writing\abc.xml	Dropdown menu
scanSpeed	4	Spin box
highlightStyle	InvertKeys	Dropdown menu
settingsFile	settings\internal_tcp.kst	Dropdown menu
windowDecorated	<input checked="" type="checkbox"/>	Checkbox
resizeEnabled	<input type="checkbox"/>	Checkbox

# Project 3: a possible solution ...



# Project 4 – Tips


---

## Requirements:

- Abotic Door Opener
  - <http://abotic.com/en/>
  - DigitalOut module
  - Plug In Abotic Door Remote to **DigitalOut port 1**
- Voice recognition
  - Microphone (Check **signal level!!**)
  - **Microsoft Speech Platform Server version 11**, the **SR language** and the **TTS language pack** must be installed too.
  - <http://www.microsoft.com/en-us/download/details.aspx?id=27225>
- FS20Sender device
  - <http://www.elv.at/fs20-funkschaltsystem.html>




# Project 4 – Tips



Check the signal level of the microphone!


Select the recognition **language**. In the **mode** „always active“ the recognition is performed continuously.

**command1-commandN**: Enter voice commands




# Project 4 – Tips

- The **DigitalOut** plugin expects action strings: "set", "clear", "toggle" and "press". The command has to be followed by a comma and the port number, for example:
  - "@DIGITALOUT:set,1" or
  - "@DIGITALOUT:toggle,2".
  - The "press"-command toggles the given output port two times with a delay of 500ms.
- **StringDispatcher** can be used to construct the string.



# Project 4 – Tips

- Enter the **housecode** and **address** of the target device in the **FS20Sender** properties.
- Our FS20Receiver Power Switch is already configured to housecode 11111111 and address 1111



# Project 4 – Tips


- In the event editor of **FS20Sender** map recognized commands with FS20 actions:
  - **off**: Switch off device
  - **onLevel1-onLevel16**: Switch on device to a certain level
  - **toggle**: Toggles switching (on/off)

Events (Ctrl-E) ▾ ↕

FS20Sender.1	SpeechProcessor.1
off	recognizedCommand4 ▾
off	--- ▾
onLevel1	--- ▾
onLevel2	--- ▾
onLevel3	--- ▾
onLevel4	--- ▾
onLevel5	--- ▾
onLevel6	--- ▾
onLevel7	--- ▾
onLevel8	--- ▾
onLevel9	--- ▾
onLevel10	--- ▾
onLevel11	--- ▾
onLevel12	--- ▾
onLevel13	--- ▾
onLevel14	--- ▾
onLevel15	--- ▾
onLevel16	recognizedCommand3 ▾
onLevel16	--- ▾
onOldLevel	--- ▾
toggle	--- ▾
dimUp	--- ▾





# Project 4: a possible solution ...



# Project 5 – Tips



## Requirements:

- Arduino
  - Connect **push button** as shown in Arduino slides, use **internal pullup resistor**
  
- IRTrans device
  - Install IRTrans driver and SW from <http://www.irtrans.de/de/download/windows.php>
  - Copy “**asterics.rem**” remote configuration file to **C:\ProgramData\IRTrans\remotes**
  - Start **IRTransServer.exe**
  - Right click on IRTrans symbol in system tray and choose “**Startup Parameter**”
  - Select “**USB**” and “**Save & Restart**”




# Project 5 – Tips

- Test sending IR codes
  - Right click on IRTrans symbol in system tray
  - Select **“Send IR Code”**
  - Select **“asterics”** as remote configuration
  - Select a command
  - Click on **“Send”**



# Project 5 – Tips

- The **CellBoard** plugin can be used to make a „button-matrix“ with scanning selection.
- scanMode „**row-column**“ highlights rows first, then columns
- Selections are done via the **scanSelect** event listener.
- Enter caption and/or image for every used button
- Use a **Timer** plugin to trigger **scanMove** events periodically.

A screenshot of the 'Properties (Ctrl-P)' window in a software development environment. The window title is 'Properties (Ctrl-P)'. The component name is 'Component: CellBoard.1'. There is a search bar with the text 'Search'. The window is divided into two main sections: 'Internal Properties' and 'Properties'.  
**Internal Properties:**


- Component Name: CellBoard.1
- Component Type: asterics.CellBoard
- Component Descriptio: GUI component
- Component Class: sensor

**Properties:**

- caption: Cell Board
- rows: 4
- columns: 2
- textColor: black
- backgroundColor: darkgray
- scanColor: red
- scanMode: row-column (highlighted with a red oval)
- hoverTime: 1000
- cellText1: start
- cellImage1: ..\OSKA\images\homecontrol\next.png
- actionText1:
- cellText2: stop
- cellImage2:
- actionText2:
- cellText3: previous
- cellImage3:
- actionText3:
- cellText4: next
- cellImage4:
- actionText4:
- cellText5: vol up
- cellImage5:
- actionText5:
- cellText6: vol down
- cellImage6:
- actionText6:
- cellText7: on / off

# Project 5 – Tips

- Use default **hostname** and **port** to connect to the local IRTrans
- **prestring** selects the remote control configuration to use, note the **comma** at the end !!
- **send1-sendN:**  
Enter the configured / stored IR codes


A screenshot of the 'Properties (Ctrl-P)' window in a software development environment. The window title is 'Properties (Ctrl-P)'. The main content area shows 'Component: IrTrans.1'. Below this, there is a search bar and two expandable sections: 'Internal Properties' and 'Properties'. The 'Internal Properties' section is expanded and shows a table with the following data:

Component Name	IrTrans.1
Component Type	asterics.IrTrans
Component Description	Infrared Transmission Module, interfaces to the
Component Class	actuator

The 'Properties' section is also expanded and shows a table with the following data:

hostname	localhost
port	21000
prestring	snd asterics,
send1	rc_play
send2	rc_stop
send3	rc_back
send4	rc_forward
send5	rc_volume_up
send6	rc_volume_down
send7	rc_operate
send8	
send9	
...	...


# Project 5: a possible solution ...



# Project 6 – Tips

---

- Add the **Ponggame** actuator and select the mode „**position**“
  - Input values of 0-300 are accepted for the paddles
  - The output of the used sensors must be adjusted to match this range of 0 – 300 to control the paddles
  - A „start“ event must be connected; this event could be generated via a GUI element ( eg. **ButtonGrid** plugin)
- Add an **Accelerometer** sensor for paddle1
  - Use a **SignalTranslation** processor to map the range of the yAcc output (-8192 – 8192) to (0 – 300)


A screenshot of a software interface showing the properties of a SignalTranslation processor. The window has a title bar with a small triangle icon and the text "Properties". Below the title bar is a table with four rows and two columns. The first column contains the property names: inMin, inMax, outMin, and outMax. The second column contains the corresponding values: -8192, 8192, 0, and 300.

Properties	
inMin	-8192
inMax	8192
outMin	0
outMax	300

# Project 6 – Tips

---

- Add the **FacetrackerLK** sensor plugin
  - The output represents head movements (relative changes of position)
  - Use the nose-Y output port (for up/down movement)
  - Use a **MathEvaluator** to invert/scale the movement values as desired
- Add an **Integration** processor to accumulate the relative changes of the head movements
  - This creates an absolute position for the paddle 2 control
  - Set upperLimit and lowerLimit of the Integrator to 300/0

A screenshot of a software interface showing the 'Properties' window for an Integrator block. The window has a title bar with a small triangle icon and the text 'Properties'. Below the title bar is a table with four rows, each representing a property and its value.

Properties	
resetValue	150
upperLimit	300
lowerLimit	0
wrapAround	<input type="checkbox"/>



# Project 6: a possible solution ...

