

# AsTeRICS Hands On! Camera Mouse Solution



# The Basics

---





- AsTeRICS can be downloaded from the official site:  
<http://www.asterics.org> or  
<https://github.com/asterics/AsTeRICS/releases>
- The latest source code is available at GitHub  
<https://github.com/asterics/AsTeRICS>
- ACS is the graphical editor, which sends models with connected plugins to the ARE (the runtime environment)

# AsTeRICS approach: Graphical Model Design

## AsTeRICS Configuration Suite (ACS)

- Connect sensor, processor and actuator modules with signal paths
- Define properties and event conditions



absolutePosition	<input type="checkbox"/>
xMin	0
xMax	1366
yMin	0
yMax	768

## AsTeRICS Runtime Environment (ARE)

- Run model configuration
- Interface sensors and actuators
- Process data




- Documentation is available:
  - Use **plugin search bar** in „Components“ tab
  - Read **tooltips** displayed when hovering over graphical elements in ACS (plugin, ports, channels, properties,...)
  - Use **F1** in ACS when plugin is selected to display plugin-help !
  - User Manual with step-by-step model creation guide
  - Developer manual for creating new plugins

Let's start ACS and ARE and have a look !

# Example 1: SignalGenerator, Oscilloscope


- How to generate signals and display them
- How to trigger event above threshold value



# Example 1: SignalGenerator, Oscilloscope

---

- Create a new model in the ACS (System → **New Model**)
- Insert the „**SignalGenerator**“ sensor component  
(Components → Sensors → Simulation → SignalGenerator)
- Insert the „**Oscilloscope**“ actuator  
(Components → Actuators → Graphical User Interface → Oscilloscope)
- Connect port „out“ to „in“



# Example 1: SignalGenerator, Oscilloscope

## SignalGenerator properties

- Change frequency to **0.2 Hz**
- Also consider sendInterval, (sample rate [ms]), amplitude, phaseShift or offset

Properties	
sendInterval	20
waveForm	sine
frequency	0.2
amplitude	100
phaseShift	0
offset	0

## Oscilloscope properties


- Use default properties

Properties	
displayBuffer	3
drawingMode	autoupdate min and max
displayMode	redraw on incoming sample
drawingInterval	100
min	-100
max	100
gridColor	black
channelColor	red
backgroundColor	white
fontSize	14
caption	oscilloscope
displayGUI	<input checked="" type="checkbox"/>

# Example 1: SignalGenerator, Oscilloscope

---

- Insert the „**Threshold**“ component  
(Components → Processors → Basic Math → Threshold)
- Insert the „**EventVisualizer**“ actuator  
(Components → Actuators → Graphical User Interface → EventVisualizer)
- Connect port „out“ to „in“ of Threshold





# Example 1: SignalGenerator, Oscilloscope

## Create Event Condition

- Connect event trigger of Threshold plugin to event listener of EventVisualizer
- Assign “eventPosEdge” to “eventDisplay”

Events (Ctrl-E)	
EventVisualizer.1	Threshold.1
eventDisplay	eventPosEdge
eventDisplay	---

## Threshold properties

- Set thresholdHigh and thresholdLow to **99**

Properties	
thresholdHigh	99
thresholdLow	99
outputHigh	1
outputLow	0
operationMode	binary
eventCondition	below->above


# Example 1: SignalGenerator, Oscilloscope

- Insert the „**WavefilePlayer**“ component (Components → Actuators → Audio and Voice)
- Connect event trigger of Threshold plugin to event listener of WavefilePlayer
- Assign “eventPosEdge” to “Start”

Events (Ctrl-E)	
WavefilePlayer.1	Threshold.1
Start	eventPosEdge
Start	---

- Select wave file to play in WavefilePlayer properties


Properties	
filename	data\sounds\applause.wav



## Example 2: Mouse-Control by Head Movements

- In this example we will show how you can control the mouse cursor via head movements
- We use a webcam and the FacetrackerLK sensor plugin (which tracks face movements)
- Mouse clicks are generated via a period of inactivity (no movements)

Note: you find a similar description in the [AsteRICS User manual](#)



just  
another  
example



## Example 2: Mouse-Control - Cursor Movement

---

### How to provide headtracking-controlled mouse alternatives for computer input

- The x- and y-position of the local mouse will be controlled by the user's head movement

#### Requirements:

- A webcam

#### Remarks:


- Camera position directly in front of the user, distance 70-120cm
- No other persons face should be in the field-of-view of the camera



## Example 2: Mouse-Control - Cursor Movement


---

- Create a new model in the ACS (System → **New Model**)
- Insert the „**FacetrackerLK**“ sensor component  
(Components → Sensors → ComputerVision → FacetrackerLK)
- Insert the „**Mouse**“ actuator  
(Components → Actuators → Input Device Emulation → Mouse)
- Connect ports „noseX“ to „mouseX“ and „noseY“ to „mouseY“




## Example 2: Mouse-Control - Cursor Movement

- Mouse plugin properties:
  - adjust **xMax** and **yMax** to the desired values, e.g. screen resolution
  - enter 0 for xMax and yMax to enable automatic detection of screen resolution
  - **deselect** „absolutePosition“  
this defines that the X and Y input values are relative changes  
(The FacetrackerLK plugin only tracks relative movements)



## Example 2: Mouse-Control - Cursor Movement


- Mouse plugin, Input Port tab:



- Select „**Synchronize Inputs**“ option in the Input Port tabs for both inputs (mouseX and mouseY)


→ this will wait for both input coordinates to arrive before the mouse position is updated

- You can try out your model by uploading it to the ARE and clicking „Start Model“



## Example 2: Mouse-Control – Add Dwell Click


- No clicking function implemented yet
- → Further development: Add **dwell clicking**
- Adding processing component „**Deadzone**“  
(Components – Processors – Signal Shaping – Deadzone)
- Connect noseX/noseY outputs of the FacetrackerLK to the inX/inY inputs of the Deadzone






## Example 2: Mouse-Control – Add Dwell Click

- Use the **Deadzone** component to define a desired movement level to start or stop the timing for the dwell click
- Deadzone component fade out x/y signal values in an adjustable range and generate event trigger if the x/y values are in- or outside this range
- Parameter „**radius**“ defines this range  
→ here it is the amount of nose movement
- Leave the radius at the default value of 10




## Example 2: Mouse-Control – Add Dwell Click

- Default value 10 – movement range is set to 10 pixels from previous to current position
- Select „**Synchronize Inputs**“ option in the Input Port Riders of the Deadzone plugin for inX and inY



## Example 2: Mouse-Control – Add Dwell Click


- How can we measure a certain timespan of low movement (to generate a dwell click) ?
- → Insert a **Timer** sensor component (Sensors – Simulation – Timer)
  - this component measures time, generates events if a time period has passed, performs timing loops
- Set time-period to 1000 ms in the component's properties



## Example 2: Mouse-Control – Add Dwell Click

- Connect event trigger port of the Deadzone component (purple) to the event listener of the Timer component (green)
- Click on the event channel (connection line is now purple dotted)

→ Events can be selected from dropdown menu




## Example 2: Mouse-Control – Add Dwell Click

- Select „**enterZone**“ event for the „start“ function
- Select „**exitZone**“ event for the „stop“ and „reset“ function
  - These event connections control the Timer components
  - **If nose movements stay below selected level of 10 pixels, the Timer is started**
  - **Else, the Timer is resetted** to 0 and stopped
  - **If the movement stays low** for the full time period, the timer will generate its „**periodFinished**“ event.

Events (Ctrl-E)	
Timer.1	Deadzone.1
start	enterZone
start	---
stop	exitZone
stop	---
reset	exitZone
reset	---


## Example 2: Mouse-Control – Add Dwell Click

- Draw a channel from the Timer's event trigger port (purple) to the event listener port of the Mouse (green)



## Example 2: Mouse-Control – Add Dwell Click


- By clicking on the new channel (line is now purple dotted), the Events can be adjusted
- Assign the „periodFinished“ event to the „leftClick“ function
- After these settings, the model is usable and provides left-click

The screenshot shows a software interface window titled "Events (Ctrl-E)". It contains two columns: "Mouse.1" and "Timer.1". The "Timer.1" column has a red oval around the "periodFinished" event, which is currently selected. The "Mouse.1" column has "leftClick" selected. The rest of the events in the "Timer.1" column are represented by "---" and a dropdown arrow.

Mouse.1	Timer.1
leftClick	periodFinished
leftClick	---
middleClick	---
rightClick	---
doubleClick	---
dragPress	---
dragRelease	---
wheelUp	---
wheelDown	---
center	---
activate	---
deactivate	---
toggle	---

## Example 2: Mouse-Control – Add GUI


**Include GUI elements for adjustable mouse acceleration and dwell timing as well as different click-actions**





## Example 2: Mouse-Control – Add GUI

- Add a **Slider** component (Components – Sensors – Graphical User Interface – Slider)
- Slider properties:
  - range of value can be defined (we can leave it at 0-100)
  - set slider’s component name to “Mouse Speed”
  - set minorTickSpacing to “10”

A screenshot of a software development environment's 'Properties (Ctrl-P)' window. The window title is 'Properties (Ctrl-P)'. The component being edited is 'Component: Mouse Speed'. Below the title bar, there are icons for a list, a search, and a zoom, followed by a search input field. The properties are organized into two sections: 'Internal Properties' and 'Properties'.

Internal Properties	
Component Name	Mouse Speed
Component Type	asterics.Slider
Component Descriptio	a slider GUI element which provide
Component Class	sensor

Properties	
min	0
max	100
default	50
autosend	<input checked="" type="checkbox"/>
caption	Mouse Speed
majorTickSpacing	20
minorTickSpacing	10
alignment	horizontal
fontSize	14

## Example 2: Mouse-Control – Add GUI


---

- To modify the x/y mouse speed with the slider's values, a **MathEvaluator** processing component is needed (Components - Processors – Basic Math – MathEvaluator)
- First, the x-signal is modified:
  - Delete port connection from noseX to mouseX
  - Draw new port connection from “value” (Slider) to “inA” (MathEvaluator)
  - Draw new connection from noseX (Facetracker) to “inB” (MathEvaluator)
  - Draw a new connection from output port (MathEvaluator) to mouseX input port

## Example 2: Mouse-Control – Add GUI

- MathEvaluator properties:
  - Adjust “**expression**” property of the MathEvaluator – this defines what will be done with the inputs – in our case we will multiply inA and inB
  - Slider **position** < 50 shall **slow down mouse** speed, slider **position** > 50 shall **increase mouse** speed

→  $a/50*b$




## Example 2: Mouse-Control – Add GUI

---

- For modifying the y-signal, copy and paste (Ctrl+C & Ctrl+V) the MathEvaluator
- New connections as for the x-direction:
  - Delete port connection from noseY to mouseY
  - Draw new connection from “value” (Slider) to “inA” of the second MathEvaluator
  - Draw new connection from noseY (Facetracker) to “inB” of the second MathEvaluator
  - Draw a new connection from output port of the second MathEvaluator to mouseY input port

# Example 2: Mouse-Control – Add GUI



## Example 2: Mouse-Control – Add GUI

---


### Adding different mouse click activities via GUI by

- adding a **ButtonGrid**
  - to select next click type
- informing mouse element about the next desired mouse click
  - by sending an “action string” to the Mouse element

**Action strings** contain commands which are understood by a number of specialized actuator elements. These string contain the addressed component and the desired command  
e.g. “@MOUSE:nextclick,right”

## Example 2: Mouse-Control – Add GUI

- Add the **ButtonGrid** component (Component – Sensors – Graphical User Interface – ButtonGrid)
- ButtonGrid properties:
  - Set “buttonCaption” properties of button 1, 2 and 3 to
    - RightClick
    - DoubleClick
    - DragClick
  - Choose “horizontalOrientation” for the ButtonGrid
  - Set a desired caption, e.g. “Next click is a”

A screenshot of a software development environment's Properties window for a ButtonGrid component. The window has a title bar that says "Properties" and a small arrow icon on the left. It contains a table with two columns: the first column lists properties and the second column shows their values. The properties listed are caption, horizontalOrientation, buttonCaption1, buttonCaption2, buttonCaption3, buttonCaption4, and buttonCaption5. The caption is "Next mouse click is a", horizontalOrientation is checked with an 'X' in a box, and the first three buttonCaption properties are set to "Right Click", "Double Click", and "Drag Click" respectively. The last two buttonCaption properties are empty.

Property	Value
caption	Next mouse click is a
horizontalOrientation	<input checked="" type="checkbox"/>
buttonCaption1	Right Click
buttonCaption2	Double Click
buttonCaption3	Drag Click
buttonCaption4	
buttonCaption5	

## Example 2: Mouse-Control – Add GUI

---

- Add the **StringDispatcher** component  
(Component – Processors – Event and String Processing – StringDispatcher)
  - Translates incoming events into outgoing strings
  - If buttons are pressed, desired action strings are generated for the Mouse components
- Connect event trigger port of ButtonGrid (purple) to the event listener port of the StringDispatcher (green)
- Click on the event channel and attach
  - button1 to dispatchSlot1
  - button2 to dispatchSlot2
  - button3 to dispatchSlot3



## Example 2: Mouse-Control – Add GUI


---

- Define the strings for slot1-slot3 in the StringDispatcher properties:
  - slot1(button1): “@MOUSE:nextclick,right”
  - slot2(button2): “@MOUSE:nextclick,double”
  - slot3(button3): “@MOUSE:nextclick,drag”
- Connect output port of the StringDispatcher to the “action” input port of the mouse

## Example 2: Mouse-Control – Add GUI

---

- Change to GUI Designer tab to define a desired position for the ButtonGrid
  - Possible positions:




- "Upload Model" and "Start Model" to try it out!


## Example 2: Mouse-Control – Improve Jittering

### Further improvement tip

- If the mouse cursor is jittering a lot, use an averager (Components – Processors – Basic Math – Averager) with a BufferSize of 5 after Facetracker Outputs



# Example 2: Mouse-Control – Improve Jittering



# Some (other) useful plugins

---

- Sensors

- ComputerVision: (X)FaceTrackerLK, EyeTribe, EyeX, KinectJ4K, TuioReactivation
- Key/Mouse: KeyCapture, MouseCapture
- Generic: DigitalIn, AnalogIn, MicGPI
- GUI: ButtonGrid, Slider, Cellboard
- Bio-Signal: p2\_parser
- Inertial: Acceleration, WiiMote
- Simulation: SignalGenerator, Timer
- Smartphone Sensors (Touch): OscServer
- Filesystem: LineReader, ReadEDF

- Actuators

- Audio: MidiPlayer, WavefilePlayer, MediaPlayer
- TTS: SpeechProcessor/SyntheticVoice
- Generic: DigitalOut, serialSender
- GUI: Oscilloscope, BarDisplay, EventVisualizer, TextDisplay
- HomeControl: FS20Sender, IrTrans
- InputDeviceEmulation: Mouse, Keyboard
- Filesystem: ApplicationLauncher, modelSwitcher, LineWriter, WriteEDF
- Gaming: PongGame
- PS3 Controller Emulation (HidActuator): RemoteMouse, Remote Keyboard, RemoteJoystick

\* Needs dedicated HW or SW / driver

# Some (other) useful plugins

---

- Processors
  - Audio&Voice: **SpeechProcessor**
  - BasicMath: MathEvaluator, Threshold, Averager, ConstantDispatcher, Comparator, Differentiate / Integrate
  - SignalShaping: Deadzone, SignalTranslation
  - DSP, Filter: IIRFilter
  - Event&StringProcessing: StringDispatcher, EventDispatcher, EventFlipFlop
  - Microcontroller: **Arduino**
  - OSKA
  - SignalPathways: PathSelector, MultiSource

\* Needs dedicated HW or SW / driver